
1) Скачать и установить Docker Desktop


https://docs.docker.com/docker-for-windows/install-windows-home/ (Windows10 Home)

https://docs.docker.com/docker-for-windows/install/ (Windows10 Pro, Enterprise)





2) Скачать по ссылке https://containds.com/desktop/download/ и установить ContainDS



Installation Problems on Windows


When running the installer on Windows, you may be asked to explicitly state that you trust software supplied by Ideonate Ltd.

If you see this:




Click 'More Info' then 'Run anyway' provided it lists the Publisher as 'Ideonate Ltd'.


- 3) Скачать файл curvefitter_container-20200602-213931.containds https://biocadru.sharepoint.com/:u:/r/sites/TemporaryCNstorsge/Shared%20Documents/NLRP 3/Curvefitter/curvefitter container-20200602-213931.containds?csf=1&web=1&e=XZk6wD
- 4) Запустить ContainDS и выбрать +New → File → Browse (выбрать файл curvefitter_container-20200602-213931.containds) → Start




5) Положить файл с результатами измерений в папку Workspace



6) Нажать на WEB (в браузер откроется новое окно с приложением)



7) Введите название файла с результатами измерений, который вы положили в папку Workspace



7) Вы прекрасны! Надеюсь, все работает)